


INDEX

PG 4-5 CliniGut PLAIN CATGUT SUTURES Tensile strength maintained for 7 - 10 days. Absorption complete in 90 days. Used in tissues which heal rapidly and require minimum support-ligating superficial blood vessels and suturing subcutaneous fatty tissue.	PG 6-9 CliniGut CHROMIC CATGUT SUTURES Tensile strength maintained for 10 - 14 days. Absorption complete in 90 days. Suitable for gynaecological procedures, general surgery, urology, approximation of fatty tissue and sub-cutaneous layers.
PG 10-16 CliniSorb VIOLET BRAIDED PGA SUTURES Predictable absorption 85% of tensile strength at 14 days, 50% at 21 days and 25% at 28 days, absorption complete in 56 - 70 days. Indications are dental, general skin closure, gynaecological, plastic, general, orthopedic and urology.	PG 17-18 CliniSorb Q BRAIDED PGA SUTURES UNDYED Tensile strength maintained for 7 - 10 days. Absorption complete in 42 days. Indications for superficial soft tissue and approximation of skin and soft tissue.
PG 19 CliniTro POLY VINYLIDENE FLUORIDE SUTURES A new generation of synthetic, non-absorbable multifilament suture. Excellent handling and tying properties. Inert, very smooth and virtually no memory. Used in skin closure, vascular and orthopedic surgery.	PG 20-22 CliniSilk BRAIDED SILK SUTURES Braided multifilament non-absorbable suture. Silicone coated to prevent wicking, very supple, high knot security. Indicated for soft tissue approximation and ligation.
PG 23-26 CliniLon NYLON MONOFILAMENT SUTURES Synthetic monofilament, non-absorbable suture. High tensile strength, low tissue reactivity, non-capillary. Easy handling and excellent passage through tissue and secure knotting. Used for skin closure, general soft tissue approximation and ligation.	PG 27-28 CliniEster POLYESTER BRAIDED SUTURES Synthetic braided non-absorbable suture. High tensile strength. Minimal acute inflammatory reaction. Suitable for cardiac and orthopedic surgery - tendon repairs, valve replacement.
PG 29-30 CliniSolv POLY P-DIOXANONE MONOFILAMENT SUTURES Predictable absorption 70% of tensile strength at 14 days, 50% at 56 days, Absorption complete in 180 - 210 days. Indicated for skin closure and non-permanent wound support - urology, orthopedic, general and gynaecological.	PG 31-32 CliniMono Q UNDYED PGCL RAPID MONOFILAMENT SUTURES High initial tensile strength diminishing after 14 days, 60% tensile strength at 7 days, 40% at 14 days and absorption complete in 90 days. Used for skin closure and non-permanent wound support, general, gynaecological, urology and plastic.
PG 33-34 CliniLene POLYPROPYLENE BRAIDED SUTURES Synthetic monofilament non-absorbable suture. Exceptionally smooth, pliable and soft. Totally inert, high tensile strength and very good knot security. Indicated for plastic surgery, digestive procedures, cardiac vascular and arteriotomy, general surgery and hernia repairs.	

Catalogue

HOW TO USE THIS CATALOGUE

The following sequence is used to code a suture:



A PLAIN GUT STRAND, 2/0 GAUGE, ATTACHED TO A 24mm CONVENTIONAL CUTTING NEEDLE

CODE 1 MATERIAL	CODE 2 GAUGE	CODE 3 NEEDLE LENGTH IN mm	CODE 4 POINT SHAPE
1 GUT, PLAIN - CliniGut	6 = 6/0	NEEDLE LENGTH	TF TAPER FINE
2 GUT, CHROMIC - CliniGut	5 = 5/0		TH TAPER HEAVY
3 POLYGLYCOLIC ACID - CliniSorb	4 = 4/0		TC TAPER CUTTING
Q3 POLYGLYCOLIC ACID - CliniSorb Q	3 = 3/0		TCH TAPER CUTTING HEAVY
4 TROFILENE - CliniTro	2 = 2/0		RC REVERSE CUTTING
5 SILK - CliniSilk	0 = 0		FRC FINE REVERSE CUTTING
6 NYLON - CliniLon	1 = 1		CC CONVENTIONAL CUTTING
7 POLYESTER FIBRE - CliniEster	22 = 2		FCC FINE CONVENTIONAL CUTTING
8 POLY (P-DIOXANONE) - CliniSolv			ST STRAIGHT TAPER
Q8 POLY (GLYCOLIDE-CO-CAPROLACTONE) - CliniMono Q			SC STRAIGHT CUTTING
9 POLYPROPYLENE - CliniLene			BP BLUNT POINT
10 LINEN - CliniLin			SG SUPER GLIDE (REVERSE CUTTING)

NEEDLE PROFILES

◀ TAPER CUTTING ● ROUND BODY ▲ CONVENTIONAL CUTTING ▼ REVERSE CUTTING

The International Colour Coding system is used for CliniSut sutures and is as depicted on the page on which the suture appears.

HEAD OFFICE & FACTORY

TEL: 041 373 9599
FAX: 041 373 9601
EMAIL: csdpe@akaciamedical.co.za

WWW.CLINISUT.COM

ADDRESS: Harrower Rd, North End,
Port Elizabeth
6001, RSA

POSTAL: PO BOX 63706,
Greenacres, 6057

Catalogue

CliniGut

PLAIN CATGUT SUTURES

USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
19mm	3/8 Rev. Cut	75cm					1319RC	1419RC		
24mm	3/8 Rev. Cut	75cm				1224RC				
26mm	3/8 Circle Rev. Cut	75cm					1326RC			
40mm	1/2 Taper Heavy	75cm		1140TH						



Catalogue

CliniGut			CHROMIC CATGUT SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
12mm	3/8 Rev. Cut	45cm								2612RC
16mm	3/8 Conv. Cut	75cm					2316CC			
16mm	3/8 Rev. Cut	75cm						2416RC		
16mm	3/8 Taper Fine	75cm						2416TF		
19mm	3/8 Rev. Cut	75cm					2319RC	2419RC	2519RC	
20mm	1/2 Taper Fine	75cm					2320TF			
24mm	3/8 Rev. Cut	75cm					2324RC			
25mm	1/2 Taper Fine	75cm					2325TF			
26mm	1/2 Rev. Cut	75cm				2226RC				
26mm	3/8 Rev. Cut	75cm					2326RC			
26mm	1/2 Taper Fine	75cm				2226TF				
30mm	3/8 Rev. Cut	75cm				2230RC	2330RC			
30mm	1/2 Taper Fine	75cm		2130TF	2030TF	2230TF	2330TF			
35mm	1/2 Taper Fine	75cm				2235TF				
35mm	1/2 Taper Heavy	75cm	22235TH							
36mm	1/2 Rev. Cut	75cm			2036RC	2236RC				
40mm	1/2 Rev. Cut	75cm	22240RC	2140RC	2040RC					
40mm	1/2 Rev. Cut	90cm				2240RC				

Catalogue

CliniGut			CHROMIC CATGUT SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
40mm	1/2 Taper Heavy	75cm	22240TH	2140TH	2040TH	2240TH				
40mm	1/2 Taper Heavy	90cm		2140 TH 90						
40mm	1/2 Blunt Point	75cm		2140BP						
45mm	1/2 Taper Fine	75cm				2245TF				
45mm	1/2 Taper Heavy	75cm	22245TH	2145TH						
48mm	1/2 Rev. Cut	75cm		2148RC						
48mm	1/2 Rev. Cut	90cm		2148RC90						
48mm	1/2 Taper Heavy	75cm		2148TH						
48mm	1/2 Taper Heavy	90cm		2148 TH 90						
50mm	1/2 Rev. Cut	75cm		2150RC						
50mm	1/2 Taper Heavy	75cm		2150TH						
50mm	1/2 Blunt Point	75cm		2150BP						
63mm	3/8 Blunt Point	75cm		2163BP						
	Tie x 1	150cm		C215						



Catalogue

CliniSorb			VIOLET BRAIDED PGA SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
12mm	3/8 Rev. Cut	45cm								3612RC
13mm	3/8 Rev. Cut	45cm						3413RC	3513RC	
13mm	3/8 Rev. Cut	45cm Undyed								3613RC-U
13mm	3/8 Rev. Cut	75cm							3513RC75	
13mm	1/2 Taper Fine	45cm							3513TF	3613TF
13mm	3/8 Taper Fine	60cm					3313TF			
16mm	3/8 Rev. Cut	45cm							3516RC	
16mm	3/8 Rev. Cut	75cm					3316RC			
16mm	3/8 Rev. Cut	45cm Undyed						3416RC-U		
17mm	1/2 Taper Fine	70cm					3317TF	3417TF		
17mm	1/2 Taper Fine	75cm							3517TF	
19mm	3/8 Rev. Cut	45cm				3219RC	3319RC	3419RC	3519RC	
19mm	3/8 Rev. Cut	45cm Undyed					3319RC-U	3419RC-U		
19mm	3/8 Rev. Cut	75cm					3319RC75	3419RC75	3519RC75	
19mm	3/8 Rev. Cut	75cm Undyed					3319RC75-U		3519RC75-U	
20mm	1/2 Taper Fine	75cm					3320TF	3420TF		

Catalogue

CliniSorb			VIOLET BRAIDED PGA SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
20mm	1/2 Taper Fine	75cm Undyed							3520TF-U	
22mm	1/2 Taper Fine	70cm				3222TF	3322TF			
22mm	1/2 Taper Fine	75cm					3322TF75			
24mm	3/8 Rev. Cut	70cm				3224RC	3324RC	3424RC		
24mm	3/8 Rev. Cut	70cm Undyed				3224RC-U	3324RC-U			
26mm	3/8 Conv. Cut	75cm					3326CC			
26mm	3/8 Rev. Cut	45cm						3426RC		
26mm	3/8 Rev. Cut	75cm				3226RC				
26mm	1/2 Taper Fine	70cm				3226TF	3326TF			
26mm	1/2 Taper Fine	75cm			3026TF		3326TF75	3426TF		
30mm	3/8 Rev. Cut	70cm				3230RC				
30mm	1/2 Taper Fine	70cm					3330TF			
30mm	1/2 Taper Fine	75cm				3230TF	3330TF75			
30mm	1/2 Taper Heavy	75cm			3030TH					
36mm	1/2 Rev. Cut	90cm		3136RC	3036RC		3236RC			
36mm	1/2 Taper Cut	90cm		3136TC	3036TC					

Catalogue

CliniSorb										
VIOLET BRAIDED PGA SUTURES										
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
36mm	1/2 Taper Fine	70cm			3036TF		3336TF			
36mm	1/2 Taper Fine	75cm				3236TF				
36mm	1/2 Taper Heavy	70cm		3136TH						
36mm	1/2 Taper Heavy	90cm			3036 TH 90					
40mm	1/2 Blunt Point	90cm		3140BP	3040BP					
40mm	1/2 Rev. Cut	90cm	32240RC	3140RC	3040RC	3240RC				
40mm	1/2 Taper Heavy	90cm	32240TH	3140TH	3040TH	3240TH				
45mm	1/2 Taper Heavy	75cm		3145TH						
45mm	1/2 Taper Heavy	90cm	32245 TH 90	3145 TH 90						
48mm	1/2 Blunt Point	90cm		3148BP						
48mm	1/2 Rev. Cut	75cm			3048RC					
48mm	1/2 Rev. Cut	90cm	32248RC	3148RC						
48mm	1/2 Taper Heavy	90cm		3148TH						
50mm	1/2 Rev. Cut	75cm	32250RC							
50mm	1/2 Rev. Cut	90cm								
50mm	1/2 Taper Heavy	70cm	32250TH							
50mm	1/2 Taper Heavy	90cm	32250 TH 90							

Catalogue

CliniSorb			VIOLET BRAIDED PGA SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
60mm	Straight Cut	70cm				3260SC	3360SC			
60mm	Straight Cut	70cm Undyed				3260SC-U	3360SC-U			
	Tie x 1	150cm				G3025	G3024			
	Tie x 1	135cm		G317	G316	G315	G314			
	Tie x 5	70cm				G3215				



Catalogue

CliniSorb Q			PGA RAPID BRAIDED SUTURES UNDYED							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
13mm	3/8 Rev. Cut	45cm							Q3513RC-U	
19mm	3/8 Rev. Cut	45cm					Q3319RC-U	Q3419FRC-U		
26mm	3/8 Rev. Cut	75cm					Q3326RC-U			
36mm	1/2 Rev. Cut	90cm				Q3236RC-U				
40mm	1/2 Rev. Cut	90cm		Q3140RC-U						
60mm	Straight Cutting	70cm				Q3260SC-U				



Catalogue

CliniSilk										
USP Gauge			BRAIDED SILK SUTURES							
			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
12mm	3/8 Rev. Cut	45cm								5612RC
13mm	3/8 Rev. Cut	45cm						5413RC		
13mm	3/8 Rev. Cut x 2	45cm						5413RC2		
16mm	3/8 Conv. Cut	75cm					5316CC	5416CC		
19mm	3/8 Rev. Cut	45cm				5219RC		5419RC	5519RC	
19mm	3/8 Rev. Cut	75cm					5319RC75			
20mm	1/2 Taper Fine	75cm					5320TF	5420TF		
24mm	3/8 Rev. Cut	45cm					5324RC			
26mm	3/8 Rev. Cut	45cm				5226RC				
26mm	3/8 Rev. Cut	75cm				5226RC75	5326RC75			
26mm	1/2 Taper Fine	75cm				5226TF	5326TF			
30mm	3/8 Rev. Cut	45cm		5130RC	5030RC					
30mm	3/8 Rev. Cut	75cm			5030RC75					
30mm	1/2 Taper Fine	75cm			5030TF					
34mm	1/2 Taper Fine	75cm		5134TF						
35mm	3/8 Rev. Cut	75cm				5235RC				

Catalogue

CliniSilk			BRAIDED SILK SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
40mm	1/2 Rev. Cut	75cm	52240RC	5140RC	5040RC					
40mm	1/2 Taper Heavy	75cm	52240TH	5140TH						
45mm	3/8 Rev. Cut	75cm			5045RC					
60mm	Straight Cut	75cm				5260SC				
80mm	3/8 Conv. Cut	100cm				5280CC				
90mm	3/8 Conv. Cut	100cm	52290CC	5190CC	5090CC	5290CC				
	Tie x 10	75cm				S523				
	Tie x 1	180cm		S595		S593	S592			



Catalogue

CliniLon			NYLON MONOFILAMENT SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
12mm	3/8 Rev. Cut	45cm							6512RC	6612RC
13mm	3/8 Rev. Cut	45cm						6413RC		
16mm	3/8 Conv. Cut	45cm							6516CC	
16mm	3/8 Rev. Cut	45cm						6416RC	6516RC	6616RC
16mm	3/8 Rev. Cut	75cm						6416RC75	6516RC75	
18mm	3/8 Conv. Cut	45cm							6518CC	
19mm	3/8 Rev. Cut	45cm					6319RC	6419RC	6519RC	
19mm	3/8 Rev. Cut	75cm				6219FRC	6319RC75	6419RC75		
20mm	3/8 Conv. Cut	45cm					6320CC	6420CC		
22mm	1/2 Taper Fine	75cm				6222TF				
24mm	3/8 Rev. Cut	45cm					6324RC	6424RC		
24mm	3/8 Rev. Cut	75cm					6324RC75	6424RC75		
25mm	3/8 Conv. Cut	45cm		6225CC			6325FCC			
26mm	3/8 Rev. Cut	45cm		6226RC			6326RC			
26mm	3/8 Rev. Cut	75cm		6226RC75			6326RC75			
30mm	3/8 Rev. Cut	75cm			6030RC	6230RC	6330RC			

Catalogue

CliniLon			NYLON MONOFILAMENT SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
30mm	1/2 Taper Heavy	100cm		6130TH						
35mm	3/8 Rev. Cut	75cm				6235RC				
36mm	1/2 Rev. Cut	75cm				6236RC				
37mm	1/2 Rev. Cut	75cm		6137RC						
39mm	3/8 Rev. Cut	75cm					6339RC			
40mm	3/8 Rev. Cut	75cm				6240RC				
40mm	1/2 Rev. Cut	100cm	62240RC	6140RC	6040RC					
40mm	1/2 Taper Heavy	100cm		6140TH						
40mm	1/2 Taper Heavy	240cm Loop		6140 ^{THL}						
45mm	3/8 Rev. Cut	100cm				6245RC				
45mm	1/2 Taper Cut Heavy	240cm Loop		6145TCL	6045TCL					
50mm	1/2 Rev. Cut	100cm		6150RC						
50mm	1/2 Taper Heavy	100cm		6150TH						
80mm	3/8 Conv. Cut	100cm				6280CC	6380CC			
90mm	3/8 Conv. Cut	100cm	62290CC	6190CC	6090CC	6290CC				



Catalogue

CliniEster

POLYESTER BRAIDED SUTURES

USP Gauge			5mm Tape	2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length									
17mm	1/2 Taper Fine	75cm						7317TF			
17mm	1/2 Taper Fine x 2	90cm				7217TF2					
26mm	1/2 Taper Fine	75cm				7226TF					
26mm	1/2 Taper Fine x 2	90cm				7226TF2					
37mm	1/2 Taper Heavy	75cm			7137TH						
40mm	1/2 Taper Cut	75cm		72240TC							
40mm	1/2 Taper Heavy	90cm		72240TH							
48mm	1/2 Taper Heavy x 2	50cm	CS48								



Catalogue

CliniSolv			POLY P-DIOXANONE MONOFILAMENT SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
13mm	1/2 Taper Fine	70cm								8613TF
13mm	1/2 Taper Fine X 2	75cm								8613TF2
17mm	1/2 Taper Fine	70cm					8317TF	8417TF	8517TF	
19mm	3/8 Rev. Cut	45cm					8319RC			
19mm	3/8 Rev. Cut	70cm						8419RC	8519RC	
20mm	1/2 Taper Fine	70cm						8420TF		
22mm	1/2 Taper Fine	70cm					8322TF	8422TF		
24mm	3/8 Rev. Cut	70cm				8224RC	8324RC	8424RC		
26mm	3/8 Rev. Cut	70cm				8226RC				
26mm	1/2 Taper Fine	70cm				8226TF	8326TF			
30mm	3/8 Rev. Cut	70cm				8230RC				
30mm	1/2 Taper Fine	70cm				8230TF	8330TF			
36mm	1/2 Rev. Cut	70cm		8136RC	8036RC	8236RC				
36mm	1/2 Taper Fine	70cm				8236TF	8336TF			
36mm	1/2 Taper Heavy	70cm			8036TH					
40mm	1/2 Rev. Cut	90cm		8140RC	8040RC					
40mm	1/2 Taper Heavy	90cm		8140TH	8040TH	8240TH				
40mm	1/2 Taper Heavy Loop	150cm		8140THL	8040THL					
48mm	1/2 Taper Heavy	90cm		8148TH						

Catalogue

CliniSolv			POLY P-DIOXANONE MONOFILAMENT SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
48mm	1/2 Taper Cut Heavy Loop	150cm		8148TCHL						
65mm	1/2 Taper Cut Heavy Loop			8165TCHL						



Catalogue

CliniMono Q

PGCL RAPID MONOFILAMENT SUTURES

USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
11mm	3/8 Rev. Cut Violet	45cm								Q8611RC
11mm	3/8 Rev. Cut Undyed	45cm								Q8611RC-U
13mm	3/8 Rev. Cut Undyed	45cm						Q8413FRC-U	Q8513RC-U	
16mm	3/8 Rev. Cut Undyed	45cm						Q8416RC-U	Q8516RC-U	
16mm	3/8 Rev. Cut Violet	45cm							Q8516RC	
19mm	3/8 Rev. Cut Violet	70cm						Q8419RC	Q8519RC	
19mm	3/8 Rev. Cut Undyed	70cm					Q8319FRC-U	Q8419RC-U	Q8519RC-U	
24mm	3/8 Rev. Cut Violet	70cm					Q8324RC	Q8424RC		
24mm	3/8 Rev. Cut Undyed	70cm					Q8324RC-U	Q8424RC-U		
26mm	3/8 Rev. Cut Violet	70cm				Q8226RC				
26mm	3/8 Rev. Cut Undyed	70cm				Q8226RC-U	Q8326RC-U			
26mm	1/2 Taper Fine Violet	70cm					Q8326TF			


Catalogue

CliniMono Q			PGCL RAPID MONOFILAMENT SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
30mm	1/2 Taper Fine Undyed	70cm					Q8330TF-U			
36mm	1/2 Rev. Cut Violet	70cm				Q8236RC				
36mm	1/2 Rev. Cut Undyed	70cm								
40mm	1/2 Taper Heavy Violet	70cm		Q8140TH						
60mm	Straight Cutting Undyed	70cm					Q8360SC-U			
										

Catalogue

CliniLene			POLYPROPYLENE BRAIDED SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
10mm	3/8 Rev. Cut	75cm								9610RC75
13mm	3/8 Rev. Cut	75cm							9513RC75	9613RC75
13mm	1/2 Taper Fine	75cm								9613TF
13mm	3/8 Taper Fine x 2	75cm								9613TF2
13mm	1/2 Taper Fine x 2	90cm								9613TF212
16mm	3/8 Rev. Cut	75cm					9316RC			
17mm	1/2 Taper Cut	90cm							9517TC	
17mm	1/2 Taper Fine	75cm						9417TF		
17mm	1/2 Taper Fine x 2	90cm						9417TF2		
19mm	3/8 Rev. Cut	75cm						9419RC75	9519RC75	
26mm	3/8 Rev. Cut	75cm				9226RC	9326RC75			
26mm	1/2 Taper Fine	75cm			9026TF	9226TF	9326TF	9426TF		
26mm	1/2 Taper Fine x 2	90cm				9226TF2	9326TF2	9426TF2		

Catalogue

CliniLene			POLYPROPYLENE BRAIDED SUTURES							
USP Gauge			2	1	0	2/0	3/0	4/0	5/0	6/0
Needle Size	Needle Description	Strand Length								
30mm	1/2 Taper Heavy	75cm		9130TH						
40mm	1/2 Rev. Cut	100cm		9140RC						
40mm	1/2 Taper Heavy	100cm		9140TH	9040TH					
60mm	Straight Cutting	75cm				9260SC				
										



Catalogue

Bonewax			
CW810	Sterile Bone Wax		2.5g

Steel Wire

USP Gauge			6	5	4
Needle Size	Needle Description	Strand Length			
48mm	1/2 Circle Conv. Cut x 4	45cm		CSS548CC	



Catalogue

BRANCHES:

Western Cape

Kiepersol Crescent

Atlas Gardens

Contermanskloof Road

Table View, 7441

Tel: +27 21 521 4100

Eastern Cape

25 Harrower Road

Kensington

Port Elizabeth, 6001

Tel: +27 41 373 9599

Gauteng

364B Olympic Duel Avenue

North Riding, 2194

Johannesburg

Tel: +27 10 003 7900

E-MAILS

Local

csd@akaciamedical.co.za

International

export@akaciamedical.co.za